



- Perfusion fluid for microdialysis
- Easy and safe to use

Perfusion Fluid T1 for Microdialysis

Perfusion Fluid T1 is an isotonic sterile perfusion fluid especially developed for Microdialysis. It is conveniently packed in glass ampoules, 5 mL, ready to use. Perfusion Fluid T1 is approved as a medical device according to the CE and the FDA medical directive

Ordering information

Features Isotonic, steril
Contents NaCl 147 mmol (3.38 g)
 KCl 4 mmol (156 mg)
 CaCl₂ 2.3 mmol (92mg)

Total Chloride contents:
 155.6 mmol/L

Osmolality 290 mOsm/kg
Ph ca.6

Package Glass ampoules, 10 x 5 mL
Ref. No. Perfusion fluid T1 P000034



CE marked according to the Medical Device Directive, 93/42/EEC



Sterilised by β -radiation



Storage temperature: 4-25 °C



Single use only

